


DAY2

Introduction to NGS technologies and library preparation

Chiara Batini, cb334@le.ac.uk






<http://www.ncbi.nlm.nih.gov/>
<https://gold.jgi.doe.gov/>



from <http://sulab.org/2013/06/sequenced-genomes-per-year>

2004: NHGRI launches '\$1000 genome' grants




Sequencers: how they were and what they were used for

platform	# reads	max read length	
Roche 454 – 2005	200K	110bp	<i>de-novo</i> assembly metagenomics
Illumina (Solexa) GA – 2006	30M	35bp	RNA-seq
ABI SOLiD – 2007	100M	35bp	ChIP-seq
IonTorrent PGM – 2010	3M	100bp	re-sequencing amplicon-seq




Next Generation Genomics: World Map of High-throughput Sequencers

☒ Show all platforms ☐ 454 ☐ HiSeq ☐ HiSeq X Ten ☐ Illumina GA2 ☐ Ion Torrent ☐ MiSeq ☐ MinION ☐ NextSeq ☐ PacBio ☐ Polonator ☐ Proton ☐ SOLiD ☐ Ser



Sequencers: size and output



Sequencers: read length

Ion PGM




Ion Proton



Illumina HiSeq



Illumina MiSeq



Roche 454 GS FLX+




250-300bp


400bp

1kb

Sequencers: the real-time generation



Pacific Biosciences (PacBio) - 2010



Nanopore MinION - 2014

Ion PGM Ion Proton



Illumina MiSeq Illumina HiSeq



Roche 454 GS FLX+



PacBio



Nanopore MinION


250-300bp

400bp

1kb

>10kb


Sequencers: size and output



Illumina sequencing technology (1)


Sequencing by synthesis, 4 modified dNTPs together, imaging support: slide

Figure 3: Attach DNA to Surface




Bind single-stranded fragments randomly to the inside surface of the flow cell channels.

Figure 4: Bridge Amplification



Add unlabeled nucleotides and enzyme to initiate solid-phase bridge amplification.


Figure 5: Fragments Become Double Stranded



The enzyme incorporates nucleotides to build double-stranded bridges on the solid-phase substrate.


Illumina sequencing technology (2)

Figure 6: Denature the Double-Stranded Molecules




Denaturation leaves single-stranded templates anchored to the substrate.

Figure 7: Complete Amplification



Several million dense clusters of double-stranded DNA are generated in each channel of the flow cell.


Figure 8: Determine First Base



The first sequencing cycle begins by adding four labeled reversible terminators, primers, and DNA polymerase.


Illumina sequencing technology (3)

Figure 9: Image First Base



After laser excitation, the emitted fluorescence from each cluster is captured and the first base is identified.

Figure 10: Determine Second Base




The next cycle repeats the incorporation of four labeled reversible terminators, primers, and DNA polymerase.

Figure 11: Image Second Chemistry Cycle



After laser excitation, the image is captured as before, and the identity of the second base is recorded.

Figure 12: Sequencing Over Multiple Chemistry Cycles





The sequencing cycles are repeated to determine the sequence of bases in a fragment, one base at a time.


IonTorrent sequencing technology (1)

Sequencing by synthesis, 1 non-modified dNTP at the time, pH detection

support: bead/chip



IonTorrent sequencing technology (2)



Sequencers: error rates and patterns

Instrument	Primary Errors	Single-pass Error Rate (%)	Final Error Rate (%)
ABI 3730xl (capillary)	substitution	0.1-1	0.1-1
Roche 454 – All models	indel	1	1
Illumina – All models	substitution	~0.1	~0.1
Ion Torrent – all chips	Indel	~1	~1
Oxford Nanopore	deletions	≥4*	4*
PacBio RS	Indel	~13	≤1

from <http://www.molecularrecologist.com/next-gen-fieldguide-2014/>


NGS applications

Genomic DNA:

- WGS (whole genome sequencing): population-scale; single-cell
- WES (whole exome sequencing)
- custom enrichment
- amplicon-seq
- RAD-seq (restriction-site-associated DNA): very useful for organisms lacking reference genomes

De novo assembly

Overlap reads by sequence, into a single contiguous sequence (Contig)




Fill gaps using known contig end sequences



Finished chromosome sequence



Re-sequencing



What does DNA-seq tell you


- Calling of known polymorphisms
- Identification of novel polymorphism (SNPs, indels etc)
- Genomic rearrangements, large deletions and insertions
- What the genome of your sample looks like (genome assembly, chromosome assembly)

NGS applications


RNA-seq:

- strand-specific RNA-seq protocols
- single-cell transcriptomics
- Fluorescent in situ RNA seq (FISSEQ): it allows the characterization of the single-cell transcriptomics and the localization of the transcripts within the cell
- CaptureSeq
- Native elongating transcript seq (NET-seq)


RNA-seq



ChIP-seq



Location-based techniques
ChIP-seq: protein-DNA
interactions; ChIP-exo
(down to nucleotide level)
protein-RNA, RNA-DNA,
DNA-DNA interactions



from <http://www.mdpi.com/2079-7737/2/1/378/htm>

Diagnostic tools in clinical laboratories; spreading in forensics


Illumina HiSeq X Ten : a human genome for less than \$1000?
 Personalized medicine is just around the corner...

So, which is the best technology?

<https://docs.google.com/spreadsheets/d/1eoWVtwKvnbV8MaRb9hzCSpGtEjazO0XTt3a8TbPG4dY/edit?pli=1#gid=0>

Good support and communication are essential

So, how do we get our DNA ready to be sequenced then?



DNA source and quantity

non-degraded, RNA-free, high quality DNA

quality can be checked with nanodrop, specific ratios are suggested:

260/280: 1.7-2.0

260/230: >2

for quantification other methods are suggested [Qubit, picogreen]

whole genome

2ug for paired end, 20ug for mate pair

Exome

6ug or more

custom enrichment

min 5ug, more depending on the target size

amplicon seq

a few ng per amplicon may be enough, since PCR will
increase the concentration

DNA source and quantity


- whole genome
- exome
- custom enrichment
- amplicon seq
- RNA


library prep

- shearing
- adapter ligation and barcoding
- size selection
- amplification
- quantification

sequencing

clonal amplification






modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

Two main methods are used:

- Sonication: hydrodynamic shearing using acoustic energy; bubbles are formed in solution, when they explode they break the DNA
- Enzymatic reaction: enzymes randomly cut the dsDNA in a time-dependent manner

library prep – adapter ligation and barcoding





modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

- adapters: 30-50bp fragments which contain primer sites for amplification and are needed to link the fragment with the support (slide, bead)
- barcodes/indexes: 6-10bp fragments which carry a unique sequence; they are used to distinguish samples run in the same lane/chip
- Adapters and barcodes are combined into one fragment


library prep – size selection

- Needed to select fragments of the right length from a mixture of fragments generated by the fragmentation step: it is usually 200-300bp for paired end sequencing with Illumina, but it may vary depending on the read length for 454 and IonTorrent



library prep – amplification


- Amplification: needed to increase the concentration of the fragments which positively incorporated adapters



modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

library prep –quantification

- Quantification: needed to tune the quantity of template for the run – usually performed with the Agilent Bioanalyzer



Overall Results for sample 3 : lib_pool3 1:15 01

Number of peaks found: 2 Corr. Area 1: 695.4
Noise: 0.2

Peak table for sample 3 : lib_pool3 1:15 01

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Area
1	35	125.00	5,411.2	Lower Marker	83.8
2	282	335.56	1,802.2		460.4
3	10,380	75.00	10.9	Upper Marker	278.6
4	12,706	0.00	0.0		6.2

DNA source and quantity

- whole genome
- exome
- custom enrichment
- amplicon seq
- RNA


library prep

- shearing
- adapter ligation and barcoding
- size selection
- amplification
- quantification

clonal amplification

- ePCR (IonTorrent/454)
- clustering (Illumina)

sequencing



clonal amplification – ePCR (Roche 454/Ion Torrent)




modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

emulsion PCR: ideally each drop contains a single fragment, PCR reagents and a bead with primers [beads may happen to be polyclonal, will be discarded after sequencing]

each bead will have one million copies of each fragment on its surface

an enrichment step is needed afterwards in order to increase the concentration of beads with attached fragments


clonal amplification – clustering Illumina



modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

isothermal amplification: temperature is kept constant but reagent are cycled to perform amplification, not very efficient (on purpose)

clonal amplification – clustering Illumina



modified from Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

clusters usually contain around 1000 identical copies of a single template

optimal cluster density will provide optimal sequencing output

A few videos...

Illumina

<https://www.youtube.com/watch?v=womKfikWlxM>

IonTorrent

<https://www.youtube.com/watch?v=ZL7DXFPz8rU>

PacBio

<https://www.youtube.com/watch?v=NHCI8PtYCFc>

Oxford Nanopore

<https://www.youtube.com/watch?v=3UHw22hBpAk>

useful references and links

Mardis 2008 Annu Rev Genomics Hum Genet 9:387-402

Metzker 2010 Nat Rev Genet 11:31-46

<http://seqanswers.com/forums/index.php>



<http://www.frontlinegenomics.com/1649/next-generation-sequencing-how-and-why-we-got-here/>

<http://thewestheimerinstitute.org/pubs/The%20challenges%20of%20sequencing%20by%20synthesis.pdf>

<http://www.molecular ecologist.com/next-gen-fieldguide-2014/>

Thanks to:

Pierpaolo Maisano Delser

Rita Neumann

James Eales



Local
Realignment

Base quality
recalibration

Mapping

Duplicate
removal

Variant
Calling

Data QC

```
graph TD; A([Data QC]) --> B([Mapping]); B --> C([Local Realignment]); C --> D([Base quality recalibration]); D --> E([Duplicate removal]); E --> F([Variant Calling]);
```

Data QC

Mapping

Local
Realignment

Base quality
recalibration

Duplicate
removal

Variant
Calling